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## BIOLOGICAL SUBSTANCE SPILL RESPONSE KIT USAGE INFORMATION

1- Personal protective equipment is worn.

2- Viruses or biological wastes are neutralized by spraying surface disinfectant on biological spillage.

3- By placing an absorbent pad as much as the waste amount on the spillage, the waste is absorbed thoroughly.

4- Remains in the area are collected with the help of dustpan and scraper by pouring the powder absorbant and solidifying.

5- Solidified residual wastes, used pads, saucages are put into the hazardous waste bag.

6- Cleaning is terminated by cleaning the area thoroughly with a cleaning cloth and sponge..

7- After all the waste is removed and the places are cleaned, the waste area and its surroundings are disinfected with surface disinfectant.

8- Hands are thoroughly cleaned with hand and skin antiseptic and procedure is terminated.

9- All materials used, including personal protective equipment used in cleaning are placed in a hazardous waste bag. Waste bag is taken to the waste storage area by Waste Officer.  
Waste Officer Code: 4217

- In case of emergency call **1122** immediately.
- If anything you need, contact the lab administrator.
- If you have questions about spill kit, don't hesitate to mail [isg@ku.edu.tr](mailto:isg@ku.edu.tr).